


| | | |
|---|---|------------------------|
|  | <i>Document Title: 4.3.2 LC-Cannabinoid Confirmation in Blood</i> | |
| | <i>Document #:</i> | <i>Page 1 of 8</i> |
| | <i>Revision #: 8</i> | <i>Published Date:</i> |
| | <i>Document Manager: Nicholas Fillinger</i> | <i>Approved By:</i> |

4.3.2.1 Title

LC-Cannabinoid Confirmation in Blood

4.3.2.2 Purpose and Scope

The purpose of this confirmatory method is to utilize LC-MS/MS technology to provide quantitative results for Δ^9 -tetrahydrocannabinol (Delta-9 THC), Δ^8 -tetrahydrocannabinol (Delta-8 THC), 11-hydroxy- Δ^9 -tetrahydrocannabinol (11-Hydroxy Delta-9 THC), and 11-nor-9-carboxy- Δ^9 -tetrahydrocannabinol (Delta-9 Carboxy THC). For ease, and space considerations, the analytes are abbreviated $\Delta 9$ -THC, $\Delta 8$ -THC, $\Delta 9$ -THC-OH, and $\Delta 9$ -THC-COOH throughout this SOP.

4.3.2.3 Limitations

Interference studies were conducted with a combination of approximately two-hundred-fifty over the counter, prescription, and illicit drugs. While this study included many drugs that are expected to be encountered in routine casework, it did not include every possible known substance that an individual could be exposed to. The drugs included in the interference study did not cause any false positive drug results for the analytes of interest.

The presence of Δ^8 -hydroxy-tetrahydrocannabinol and Δ^8 -carboxy-tetrahydrocannabinol may interfere with Δ^9 -hydroxy-tetrahydrocannabinol and Δ^9 -carboxy-tetrahydrocannabinol, respectively. The presence of exo-tetrahydrocannabinol may interfere with Δ^9 -tetrahydrocannabinol. See 4.3.2.12 for additional details.

Carryover was assessed after analyzing standards that were ten times higher than standard six. No carryover existed in blank samples analyzed immediately afterward. To be conservative, analyte/internal standard area ratios higher than five times standard six shall result in the reinjection of the next casework sample.

4.3.2.4 Specimen Criteria

500 μ L blood

4.3.2.5 Safety

blood – mitigate biohazard risk by employing universal biohazard precautions


acetonitrile – see MSDS, mitigate chemical hazard risk by taking appropriate precautions

formic acid – see MSDS, mitigate chemical hazard risk by taking appropriate precautions

methanol – see MSDS, mitigate chemical hazard risk by taking appropriate precautions

methyl tert-Butyl Ether (MTBE) – see MSDS, mitigate chemical hazard risk by taking appropriate precautions

THIS DOCUMENT IS UNCONTROLLED IF IN HARD COPY FORM
MICHIGAN STATE POLICE FORENSIC SCIENCE DIVISION

| | | |
|---|---|------------------------|
|  | <i>Document Title: 4.3.2 LC-Cannabinoid Confirmation in Blood</i> | |
| | <i>Document #:</i> | <i>Page 2 of 8</i> |
| | <i>Revision #: 8</i> | <i>Published Date:</i> |
| | <i>Document Manager: Nicholas Fillinger</i> | <i>Approved By:</i> |

4.3.2.6 Equipment

LCMSMS-21-1
 Poroshell 120, EC-C18, 2.1x100mm, 2.7µm column
 Poroshell 120, UHPLC Guard, EC-C18, 2.1mm
 1 mL capacity Biotage ISOLUTE+ SLE cartridge
 Centrifuge
 Pipettes
 Repeater pipette
 Vortex mixer
 Solvent evaporator
 Vacuum pump
 Glass block vacuum chamber
 Extraction cartridge 12/24 manifold
 Clean through extraction tip
 16x100 mm glass culture tubes
 Autosampler vial with glass insert
 Autosampler vial cap
 Formic acid – LCMS grade
 Acetonitrile – LCMS grade
 Deionized water – LCMS grade
 MTBE

4.3.2.7 Sample Preparation and Procedure

Prepare extraction manifold by labeling a 1 mL capacity SLE cartridge and placing cartridge onto the extraction manifold equipped with disposable extraction tips. Place a collection tray with appropriately labeled 16x100 mm collection tube into the glass block vacuum chamber.

Note: ensure that each disposable extraction tip is lined up with the corresponding collection tube.

To all appropriately labeled 16x100 mm glass tubes, add 500 µL of 0.1% formic acid in water. Add 20 µL of internal standard solution to all culture tubes, including the negative control. Add 20 µL of standards 1-6, and the low, medium, and high control to appropriate culture tubes.

Prepare an interference standard by adding 20 µL of interference mix to 500 µL of blank blood. The interference mix includes:


- 1000 ng/mL CBD
- 2.5 ng/mL Δ8-THC-COOH
- 1 ng/mL Δ8-THC-OH

Prepare a matrix negative control by adding 500 µL of blank blood (no internal standard) to appropriate culture tube. A matrix negative is not necessary if the lot# of blank blood being used for calibrators and controls has previously been analyzed in this confirmatory assay. It is the responsibility of the individual performing the confirmatory analysis to verify this information.

Add 500 µL of blank blood to standard and control culture tubes.

Add 500 µL of unknown blood to culture tubes.

THIS DOCUMENT IS UNCONTROLLED IF IN HARD COPY FORM
 MICHIGAN STATE POLICE FORENSIC SCIENCE DIVISION

| | | |
|---|---|------------------------|
|  | <i>Document Title: 4.3.2 LC-Cannabinoid Confirmation in Blood</i> | |
| | <i>Document #:</i> | <i>Page 3 of 8</i> |
| | <i>Revision #: 8</i> | <i>Published Date:</i> |
| | <i>Document Manager: Nicholas Fillinger</i> | <i>Approved By:</i> |

Vortex each calibrator, control, and unknown to mix sample with diluent.
 Pipette 750 µL of each calibrator, control, and unknown into the labeled SLE cartridge.
 Apply a short pulse of vacuum (approximately 1-2 seconds) to initiate the sample loading into the SLE cartridge. Allow gravity or slight individual positive pressure to load the samples into the SLE cartridge.

Note: wait for each sample to completely descend into the SLE cartridge before moving onto the next step.

Allow cartridge sorbent 5 minutes to absorb the sample.
 Add 3 mL of MTBE to each SLE cartridge and allow 5 minutes for gravity to begin eluting the sample.
 Add a second 3 mL aliquot of MTBE to each SLE cartridge and allow 5 minutes for gravity to begin eluting the sample.

Note: after the second 5 minute elution, use minimal vacuum pressure or pulse the vacuum pump to assist with the elution of the sample. Do not exceed approximately 1 mL elution flow per minute.

Place collection tubes into solvent evaporator (<40°C) and dry with nitrogen (approximately 4-6 psi) until dry (approximately 5-10 minutes).
 Reconstitute by adding 40 µL of 0.1% formic acid in acetonitrile to each culture tube, vortex. Add 40 µL of 0.1% formic acid in water to each culture tube, vortex.
 Transfer each sample into an appropriately labeled LC autosampler vial and cap.
 Centrifuge for 10 minutes at 3500 rpm.
 Run on instrument.

Note: validation demonstrated that the analytes in this method are stable for 3 days after extraction, when stored in 15 °C. Samples not analyzed within 3 days of extraction must be reextracted.

4.3.2.8 Performance Characteristics


Standards have been validated at the following concentrations (ng/mL):

| Analyte | Standard 1 | Standard 2 | Standard 3 | Standard 4 | Standard 5 | Standard 6 |
|-------------|------------|------------|------------|------------|------------|------------|
| Δ9-THC | 1 | 5 | 30 | 50 | 70 | 100 |
| Δ8-THC | 1 | 5 | 30 | 50 | 70 | 100 |
| Δ9-THC-COOH | 5 | 10 | 25 | 50 | 75 | 100 |
| Δ9-THC-OH | 1 | 4 | 8 | 12 | 16 | 20 |

Controls have been validated at the following concentrations (ng/mL):

| Analyte | LOW Control | MEDIUM Control | HIGH Control |
|-------------|-------------|----------------|--------------|
| Δ9-THC | 3 | 40 | 80 |
| Δ8-THC | 3 | 40 | 80 |
| Δ9-THC-COOH | 8 | 40 | 80 |
| Δ9-THC-OH | 2 | 10 | 18 |

THIS DOCUMENT IS UNCONTROLLED IF IN HARD COPY FORM
 MICHIGAN STATE POLICE FORENSIC SCIENCE DIVISION

| | | |
|---|---|-----------------|
|  | Document Title: 4.3.2 LC-Cannabinoid Confirmation in Blood | |
| | Document #: | Page 4 of 8 |
| | Revision #: 8 | Published Date: |
| | Document Manager: Nicholas Fillinger | Approved By: |

Note: Instructions for calibrator and control preparation can be found in form FS-127.

Internal Standards have been validated at the following concentrations (ng/mL):

| Analyte | Concentration |
|----------------|---------------|
| Δ9-THC-D3 | 50 |
| Δ8-THC-D3 | 50 |
| Δ9-THC-COOH-D3 | 50 |
| Δ9-THC-OH-D3 | 10 |

Note: Instructions for Internal Standard preparation can be found in form FS-126.

Limits of Detection (LOD), Lower Reporting Limits (LRL), Lower Limits of Quantitation (LLOQ), and Upper Limits of Quantitation (ULOQ) have been validated at the following concentrations (ng/mL):

| Analyte | LOD | LRL | LLOQ | ULOQ |
|-------------|-----|-----|------|------|
| Δ9-THC | 0.4 | 0.5 | 1 | 100 |
| Δ8-THC | 0.4 | 0.5 | 1 | 100 |
| Δ9-THC-COOH | 2 | 2.5 | 5 | 100 |
| Δ9-THC-OH | 1 | 1 | 1 | 20 |

Note: the LOD for Δ9-THC, Δ8-THC, and Δ9-THC-COOH may be lower than what is listed in this table. As no concentrations below the LRL will be reported, validation studies did not go any lower.

Note: the LRL is determined as 50% of the area ratio between analyte and internal standard of calibrator 1. The concentrations provided in the above table are estimates of the LRL. The LRL for Δ9-THC-OH is equivalent to Standard 1.


4.3.2.9 Analytical Parameters

Instrument Operating Parameters are in TX-PM 5.15 Instrument Operating Parameters.

Instrument Maintenance

- Check to ensure that the following maintenance has been completed:
- Curtain plate has been cleaned within the last 30 days
- Guard column has been changed within the last 2 months
- Mobile phases have been prepared within the last month and are of sufficient volume to complete the analytical run
- HPLC is connected to MS
- Correct mobile phases are being used
- Correct column is in place
- Waste container is not full
- Pump valves are closed

THIS DOCUMENT IS UNCONTROLLED IF IN HARD COPY FORM
MICHIGAN STATE POLICE FORENSIC SCIENCE DIVISION

| | | |
|---|---|------------------------|
|  | <i>Document Title: 4.3.2 LC-Cannabinoid Confirmation in Blood</i> | |
| | <i>Document #:</i> | <i>Page 5 of 8</i> |
| | <i>Revision #: 8</i> | <i>Published Date:</i> |
| | <i>Document Manager: Nicholas Fillinger</i> | <i>Approved By:</i> |

Needle rinse has sufficient volume

Calibration Models

| Analyte | Curve Type | Weighting | Origin | Transition 1 | Transition 2 | Internal Standard |
|-------------|------------|---------------|--------|---------------|---------------|---------------------------------|
| Δ9-THC | Quadratic | Inverse (1/x) | Ignore | 315.1 / 193.1 | 315.1 / 123.0 | Δ9-THC-D3 318.1 / 123.0 |
| Δ8-THC | Quadratic | Inverse (1/x) | Ignore | 315.1 / 193.1 | 315.1 / 123.1 | Δ8-THC-D3 318.1 / 123.0 |
| Δ9-THC-COOH | Linear | Inverse (1/x) | Ignore | 343.0 / 299.1 | 343.0 / 191.0 | Δ9-THC-COOH-D3 346.0 / 194.0 |
| Δ9-THC-OH | Linear | Inverse (1/x) | Ignore | 331.1 / 193.1 | 331.1 / 105.1 | Δ9-THC-OH-D3 334.1 / 196.1 |

Note: The retention time order of the analytes is Δ9-THC-OH, Δ9-THC-COOH, Δ9-THC, and Δ8-THC.


4.3.2.10 Data Analysis

LCMSMS-21-1 Cannabinoid Flagging Rules and Calculated Columns

When processing a batch with Sciex OS the results table fields listed below shall comply with the following numbering format:

| Column name | Number format |
|--------------------------|---------------|
| Accuracy % | 0.00 |
| Actual Concentration | 0.00 |
| Area | 0.000e0 |
| Area Ratio | 0.0000 |
| Calculated Concentration | 0.000 |
| Exp. IR High | 0.000 |
| Exp. IR Low | 0.000 |
| Exp. RRT High | 0.000 |
| Exp. RRT Low | 0.000 |
| Expected Ion Ratio | 0.000 |
| Expected RRT | 0.000 |
| Ion Ratio | 0.000 |
| ISTD Area | 0.000e0 |
| Mean ISTD Area | 0.000e0 |

THIS DOCUMENT IS UNCONTROLLED IF IN HARD COPY FORM
MICHIGAN STATE POLICE FORENSIC SCIENCE DIVISION

| | | |
|---|---|------------------------|
|  | <i>Document Title: 4.3.2 LC-Cannabinoid Confirmation in Blood</i> | |
| | <i>Document #:</i> | <i>Page 6 of 8</i> |
| | <i>Revision #: 8</i> | <i>Published Date:</i> |
| | <i>Document Manager: Nicholas Fillinger</i> | <i>Approved By:</i> |

| | |
|-----------------------|--------|
| Relative RT | 0.000 |
| Reporting Limit | 0.0000 |
| Standard 1 Area Ratio | 0.0000 |

Acceptance Criteria of Calibrators

Calibrator values shall be within $\pm 20\%$ of the target value. The r^2 value for calibration curves shall be ≥ 0.990 .

Method Validation did not evaluate the impact of dropping calibration points on the linear curves of $\Delta 9$ -THC-OH or $\Delta 9$ -THC-COOH. No calibration points may be dropped.

Acceptance Criteria of Controls

Control values shall be within $\pm 20\%$ of the target value. All quantitative controls are required to be within acceptance criteria. The negative control shall not contain the analytes of interest at a reportable concentration.

The interference standard shall not contain the analytes that are being reported in this confirmatory assay. The peak height of $\Delta 8$ -THC-COOH shall be used as a threshold to determine the reportability of $\Delta 9$ -THC metabolites.

Acceptance Criteria of Internal Standards

Internal Standard area counts shall be within acceptance criteria of the mean of the batch, specific to each IS. The following acceptance criteria shall be used:

| Analyte | Acceptance |
|-------------------------|------------|
| $\Delta 9$ -THC-D3 | $\pm 30\%$ |
| $\Delta 8$ -THC-D3 | $\pm 40\%$ |
| $\Delta 9$ -THC-COOH-D3 | $\pm 25\%$ |
| $\Delta 9$ -THC-OH-D3 | $\pm 20\%$ |


Acceptance Criteria of Casework Samples

Chromatographic Quality Requirement

Chromatographic quality for LC/MS/MS transition data is defined as a reasonably symmetrical shaped peak consistent with those observed in calibrators and positive controls and can be differentiated from a negative control.

Retention Time Requirement

THIS DOCUMENT IS UNCONTROLLED IF IN HARD COPY FORM
MICHIGAN STATE POLICE FORENSIC SCIENCE DIVISION

| | | |
|---|---|------------------------|
|  | <i>Document Title: 4.3.2 LC-Cannabinoid Confirmation in Blood</i> | |
| | <i>Document #:</i> | <i>Page 7 of 8</i> |
| | <i>Revision #: 8</i> | <i>Published Date:</i> |
| | <i>Document Manager: Nicholas Fillinger</i> | <i>Approved By:</i> |

This assay uses a Relative Retention Time (RRT) model. The RRT is calculated by dividing the analyte retention time by the internal standard retention time. This is performed for all calibrators used in the calibration curve, averaged, and then given an acceptable value of $\pm 2.5\%$. Positive analytes must be within this range to be considered acceptable.

Ion Ratio Requirement

This assay utilizes one transition ion as a quantitative ion and a second transition ion as a confirmatory ion. The ion ratio is calculated by dividing the area count of transition 2 by the area count of transition 1. The ratio between the two ions must be within $\pm 20\%$ of the average of all the calibrators used in the calibration curve to be acceptable.

NOTE: Chromatographic and Mass Spectral Quality Control shall be evaluated for all calibrators, controls, and unknowns in an analytical batch.

4.3.2.11 Corrective Measures

If acceptance criteria are not met, consult with the unit supervisor to determine whether re-extraction or qualitative reporting will be performed. The rationale for that decision shall be recorded in the electronic case file.

4.3.2.12 Reporting

Quantitative results will be truncated and reported in ng/mL to one decimal for values <10 ng/mL. Quantitative values ≥ 10 ng/mL will be truncated and reported as a whole number.


Analytes present at a concentration $<$ (LOQ) with an area ratio of analyte to I.S. $\geq 50\%$ of Standard 1 will be reported as " $<$ (LOQ)".

Analytes present at a concentration $<$ LOQ with an area ratio of analyte to internal standard $< 50\%$ of Standard 1 will not be reported. If no analytes in this confirmatory assay are positive the report will state, "Not Detected: Cannabinoids".

If an analyte is present above the ULOQ, it will be reported as " $>$ (ULOQ)". If an analyte is present above the ULOQ and the ion ratios fail acceptance criteria, the sample may be reinjected with $1 \mu\text{L}$ to demonstrate acceptable ion ratios. The impacted analyte will remain reported as " $>$ (ULOQ)".

If $\Delta 8\text{-THC-COOH}$ is present above the interference standard threshold, $\Delta 9\text{-THC}$ metabolites cannot be reported. In those instances, the language on the reports shall state: " $\Delta 9\text{-THC-COOH}$ and $\Delta 9\text{-THC-OH}$ were not able to be reported due to an interfering substance. For additional information contact this laboratory."

Validation demonstrated that ≥ 10 ng/mL of exo-tetrahydrocannabinol interferes with $\Delta 9\text{-THC}$, at all calibration levels. When both analytes are present in the same sample, the report shall state: " $\Delta 9\text{-THC}$ was not able to be reported due to an interfering substance. For additional information contact this laboratory."

| | | |
|---|---|------------------------|
|  | <i>Document Title: 4.3.2 LC-Cannabinoid Confirmation in Blood</i> | |
| | <i>Document #:</i> | <i>Page 8 of 8</i> |
| | <i>Revision #: 8</i> | <i>Published Date:</i> |
| | <i>Document Manager: Nicholas Fillinger</i> | <i>Approved By:</i> |

4.3.2.13 References

LCMSMS-21-1 LC-Cannabinoid Final Validation Document

ANSI/ASB Standard 036

ANSI/ASB Standard 054

ANSI/ASB Standard 152

ANSI/ASB Standard 017